GUIDE CATHETER EXTENSION

Instructions For Use	2
Инструкции за употреба – Bulgarian	4
Upute za uporabu – Croatian	6
Návod k použití – Czech	8
Brugsanvisning – Danish	10
Gebruiksaanwijzing – Dutch	12
Kasutusjuhend – Estonian	14
Käyttöohjeet – Finnish	16
Mode d'emploi – French	18
Gebrauchsanweisung – German	20
Οδηγίες χρήσης – Greek	22
Használati utasítás – Hungarian	24
Istruzioni per l'uso – Italian	26
Lietošanas norādījumi – Latvian	28
Naudojimo instrukcijos – Lithuanian	30
Bruksanvisning – Norwegian	32
Instrukcja użycia – Polish	34
Instruções de utilização – Portuguese	36
Instrucţiuni de utilizare – Romanian	38
Návod na použitie – Slovak	40
Instrucciones de uso – Spanish	42
Bruksanvisning – Swedish	44
Kullanma Talimatı – Turkish	46

Guide Catheter Extension

DEVICE DESCRIPTION

The LiquiD Guide Catheter Extension is a single lumen catheter offered in a 6F and 7F size. The 150cm device has a stainless-steel shaft connected to a 15cm single lumen tube, which is used to extend traditional guide catheters. The 15cm single lumen tube catheter body contains a coil for kink resistance and radiopacity. The single lumen tube catheter body is also silicone coated for lubricity. The device has two proximal positioning marks located at 95cm and 105cm from the distal tip. The device handle is color coded to match the standard guide catheter color code.

DEVICE COMPATIBILITY

MODEL SIZE	ID	OD	COMPATIBLE GUIDE CATHETER
LiquiD 061 6F	0.061 inches (1.55mm)	0.068 inches (1.73mm)	≥6F
LiquiD 071 7F	0.071 inches (1.80mm)	0.078 inches (1.98mm)	≥7F
LiquiD 051 6F	0.051 inches (1.30mm)	0.068 inches (1.73mm)	≥6F

INTENDED USE

The LiquiD Guide Catheter Extension is intended to be used in conjunction with guide catheters to access discrete regions of the coronary and/or peripheral vasculature, and to facilitate placement of interventional devices.

CONTRAINDICATIONS

- Vessels less than 2.5mm in diameter
- Vessels in the neurovascular system and the venous system

WARNINGS

- Do not use the device if the sterile barrier has been damaged. Use of a non-sterile device may cause infection.
- For single use only. Do not reuse, reprocess, or re-sterilize the device. Reuse, reprocessing, or re-sterilization may compromise device integrity resulting in device separation or vessel injury. Reuse may also cause device contamination which may result in infection.
- Failure to abide by the warnings in this labelling may result in damage to the device coating, which may necessitate intervention or result in serious adverse events.
- Coating integrity or performance may be negatively impacted by preparation with solvents (e.g. alcohol or antiseptic).
- Use prior to the "use by" date as indicated on the label
- Never advance the LiquiD device into a vessel without a leading guidewire and without confirming distal tip location using fluoroscopic guidance. Vessel injury may result.
- Never advance the LiquiD device into a vessel with an effective diameter less than 2.5mm. Vessel injury, ischemia, and/or occlusion may result. If pressure in a vessel dampens after inserting the LiquiD device, withdraw the LiquiD device immediately.
- Never advance the LiquiD device more than 10cm beyond the guide catheter tip. The LiquiD device may become lodged in the guide catheter making it difficult to remove.

- Never advance or withdraw an intravascular device against resistance until the cause of resistance is determined by fluoroscopy. Movement of the catheter or guidewire against resistance may result in separation of the catheter or guidewire, catheter damage, or vessel injury.

PRECAUTIONS

- This device should be used only by physicians trained in percutaneous, intravascular techniques and procedures.
- Inspect the device prior to use for any bends or kinks. Use of a damaged device may result in device separation or vessel injury.
- Flush the device lumen with heparinized saline prior to use.
- Use of systemic anti-coagulation therapy should be considered to prevent or reduce clotting when any device is used in the vascular system.
- Handle the device with care during the procedure to prevent the possibility of accidental breakage, bending, or kinking.
- This device should not be torqued more than 360 degrees. Over torquing the device may result in device damage, entrapment or vessel injury.
- Manipulate the device only when observing the distal tip under fluoroscopy when in the body to reduce the risk of vessel injury.

ADVERSE EVENTS

The adverse events include, but are not limited to:

- Allergic reaction
- Arterial spasm
- Cardiac arrest
- Embolism
- Hemorrhage/hematoma
- Infection
- Ischemia
- Myocardial infarction
- Occlusion
- Stent dislodgement/damage
- Thrombus
- Vessel injury (i.e. perforation, dissection, intimal disruption)

HOW SUPPLIED

The LiquiD device is supplied STERILE using an ethylene oxide (EO) process. Do not use if the sterile barrier is damaged. Packaging is designed to maintain sterility according to labeled "use by" date. Do not use if labeling is incomplete or illegible.

ANCILLARY DEVICES

- Guide catheter
- Y-adaptor with hemostasis valve
- Guidewire
- 10cc sterile syringe for flushing
- Sterile heparinized saline

DEVICE OPERATION

Preparations for Use

1. Inspect the LiquiD device sterile barrier for damage.
2. Transfer the dispenser coil with the LiquiD device into the sterile field using standard sterile handling procedures.
3. Remove the LiquiD device from the dispenser coil.
4. Flush the LiquiD device lumen with heparinized sterile saline using a 10cc syringe via the distal tip.

Instructions for Use

1. Inspect the LiquiD device prior to use to verify that it has not been damaged during handling.

2. Backload the LiquiD device onto the guidewire and advance until the catheter is proximal of the hemostasis valve.
3. Open the hemostasis valve and advance the LiquiD device into the guide catheter.
4. Advance the LiquiD device through the guide catheter and into the desired vessel position.

WARNING: Never advance the LiquiD device into a vessel without a leading guidewire and without confirming distal tip location using fluoroscopic guidance. Vessel injury may result.

WARNING: Never advance the LiquiD device more than 10cm beyond the guide catheter tip. The LiquiD device may become lodged in the guide catheter making it difficult to remove.

WARNING: Never advance the LiquiD device into a vessel with an effective diameter less than 2.5mm. Vessel injury, ischemia, and/or occlusion may result. If pressure in a vessel dampens after inserting the LiquiD device, withdraw the LiquiD device immediately.

5. Insert the interventional device over the guidewire, through the hemostasis valve, and into the guide catheter.
6. Close the hemostasis valve. Perform the interventional procedure according to the instructions provided by the interventional device manufacturer.
7. Remove all interventional devices.
8. Remove the LiquiD device ensuring that the hemostasis valve is fully open prior to retracting the LiquiD device lumen through the hemostasis valve.
9. Dispose of device and packaging in accordance with hospital administrative and local government policy.

Notice

If a serious incident related to the LiquiD device has occurred, report the incident to the manufacturer and the local government competent authority.

	Batch code
	Catalog number
	Medical Device
	Sterilized using ethylene oxide
	EU authorized representative
	Legal manufacturer
	Prescription only
	Consult instructions for use
	Date of manufacture
	Use by date
	Do not re-use
	Keep away from rain
	Keep away from sunlight
	Non pyrogenic
	Do not use if package is damaged
	Do not resterilize
	Inner diameter
	Single sterile barrier system with outer protective packaging

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REGULATORY INFORMATION

A pdf copy of this IFU is located at (<http://www.seiglamedical.com/IFU>)

This is the LiquiD Guide Catheter Extension Devices - 6F, 7F, 6F LP

(Basic UDI-DI: 0850012180M10006F4X, 0850012180M20007F5F, 0850012180M30006FLPDV)

Summary of Safety and Clinical Performance (SSCP) location after the launch of the European Database on Medical Devices/EUDAMED (<http://ec.europa.eu/tools/eudamed>)

EC REP EU Authorized Representative

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UK RP UK Responsible Person

MDSS-UK RP LIMITED
Parkway House,
Palatine Rd, Northenden,
Wythenshawe,
Manchester M22 4DB,
United Kingdom

 Legal Manufacturer

Seigla Medical Inc.
7688 5th Street Southeast
Buffalo, MN 55313
USA
Telephone 612.615.9058

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MDSS GmbH